Owner's Manual
Manuel De L'utilisateur
Manual Del Propietario
Benutzerhandbuch

 **MAGNASSAGER**

THANE Fitness®

Contents

1	Introduction
2	Before You Begin
3	Important Safeguards
5	Product Features
7	When to Use and When Not to Use Your Maganassager
8	Special Considerations
9	How to Use the Magnassager

Introduction

Thank you for selecting the Magnassager. I hope you discover and enjoy the many benefits from the use of this product.


There are no substitutes for a professional hand massage, however, the Magnassager is the only product that simulates the actual kneading technique of massage. It was designed to provide a soothing massage without the irritating vibration that is characteristic of other massage devices. Another innovative feature is the natural permanent magnet which suspends the massage balls in place.

With the action of the rotating massage balls, the correct massage technique is similar to hand massage which incorporates a slow, but stronger pressure toward the heart and a softer return stroke away from the body. The user should allow the massage balls to freely roll over the muscles. When continuous massage is maintained, heat conducts into the deep muscle layers.

I am confident that as you use your Magnassager, you will be impressed not only with its functional design and innovative action, but you will grow to appreciate the quality components that have gone into the manufacture of the finest massage device in the industry.

Once again, thank you for selecting the Magnassager.

Sincerely,

A handwritten signature in black ink, reading "Rick E. Jones".

Before You Begin

The application of the Magnassager is unlike any massage device on the market. Most users have come to expect massage devices to vibrate or pulverize the muscle into submission. The Magnassager does not vibrate; it is intended to soothe and gently knead the muscle. Please follow the recommended techniques for best results:

- *Use slow speed settings to get used to the Magnassager's action.*
- *Maintain constant contact with the body and all three massage balls.*
- *Massage the muscle; keep clear from bony protrusions.*
- *For individuals with low muscle mass, try the Magnassager over a towel held taut.*
- *Use slow massage strokes.*
- *Let the Magnassager do the work for you.*

IMPORTANT SAFEGUARDS

FOR HOUSEHOLD USE ONLY

When using an electrical appliance, basic precautions should always be followed including the following:

READ ALL INSTRUCTIONS BEFORE USING THE "MAGNASSAGER"

DANGER

To reduce the risk of electric shock:

- Read all instructions before using the appliance.
- DO NOT use while bathing or showering.
- Do not place or store product where it can fall or be pulled into a tub or sink.
- Do not place in or drop into water or other liquid.
- Do not reach for an appliance that has fallen into water. Unplug immediately.
- Always unplug this appliance from the electrical outlet immediately after using and before cleaning.
- Do not use water to clean or sterilize appliance. For cleaning, you can use a damp cloth, **ONLY** after unplugging the appliance.

WARNING

To reduce the risk of burns, fire, electric shock or injury to persons:

- Close supervision is necessary when the Magnassager is used by, on, or near children, invalids, or disabled persons.
- Use the Magnassager only for its intended use as described in this manual. Do not use any other accessories or attachments not recommended by the manufacturer.
- Never operate the Magnassager if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or if it has been dropped into water. Return the appliance to the factory or manufacturer for examination and/or repair.
- Keep the electrical cord away from heated surfaces.

IMPORTANT SAFEGUARDS (cont.)

- Always grab by the plug to remove from outlet. DO NOT pull on cord.
- Never force the plug into an outlet.
- Never use while drowsy.
- Never block the air openings of the appliance or place it on a soft surface, such as a bed or couch where the air openings may be blocked. Keep the air openings free of lint, hair, and the like.
- Never drop or insert any objects into any opening or hose.
- Do not use outdoors or operate where aerosol (spray) products are being used or where oxygen is being administered.
- The appliance should never be left unattended when plugged in. Unplug from outlet when not in use, and before attaching or removing parts.
- Do not operate under a blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
- Do not operate near televisions, personal computers, floppy disks, compact disks, pagers, cellular phones or other sensitive electronic equipment.
- The Magnassager contains a natural permanent magnet which emits a negative pulsating magnetic field of 1000 gauss. **Individuals with pacemakers or detached defibrillators may be at risk.**
- Do not carry this appliance by the supply cord or use the cord as a handle.
- To disconnect, turn all controls to the OFF position, then remove plug from outlet.
- Never use pins, knives, tacks or any other metallic fasteners or objects near the Magnassager.
- Keep children away from this appliance. If swallowed, massage balls may become lodged in the throat.
- Do not drop or slam the Magnassager against any surface. If dropped on control knob, it may short circuit. If dropped on the massage balls, it may impact the magnet and break. Return the appliance to the factory or manufacturer for examination and/or repair.
- **Do not use the Magnassager on patients with pacemakers or internal defibrillators.**

SAVE THESE INSTRUCTIONS

Product Features

■ User Maintenance

The Magnassager comes pre-assembled. To dismantle for cleaning, make sure the speed control knob is in the OFF position and disconnect the unit from the power supply.

The stainless steel balls are suspended in place by a magnet. To remove the balls from the rotator plate, simply pull them out with your fingers or shake the unit.

After removing the stainless steel balls from the rotator plate, detach the plate from the mounting shaft by pulling down on the rotator plate.

Wash the stainless steel balls and rotator plate with an anti-bacterial soap. Do not use thinner or other solvents. Dampen a cloth with a mild soap and wipe the inside and outside of the Magnassager housing. Do not submerge Magnassager in water. Allow a few minutes for the Magnassager to dry.

Periodically, clean the air vents with a dry toothbrush to remove lint deposits.

■ Assembly

Position the rotator plate to slip on to the shaft adapter. The rotator plate should be flush with the housing after installation.

Place the stainless steel balls into the holes. The magnet should pull them into place.


■ Operation

Note: Always check for metallic objects attached to the bottom of device. The Magnassager's switch should be in the OFF position when you do this.

While holding the Magnassager, plug in wall adapter. Turn the switch clockwise to the ON position. Turn knob to desired speed. Apply the Magnassager to massage site.

■ Product and Cord Storage When Appliance Is Not in Use.

The Magnassager should be unplugged and stored after each use away from heat and water. Cord should not be allowed to hang over edge of counter or table or touch hot or wet surfaces. The cord should be coiled after each use and stored.

Product Features (cont.)

■ Cord Care Instructions

When appliance is in use check cord periodically to make sure it has not become twisted or tangled. If this occurs gently untwist cord and remove tangles or knots.

DANGER

This product is factory equipped with a specific electric cord and plug to permit connection to a proper electric circuit. Make sure that the product is connected to an outlet having the same configuration as the plug. No adapter should be used with this product. Do not modify the plug provided – if it will not fit the outlet, have the proper outlet installed by a qualified electrician. If the product must be reconnected for use on a different type of electric circuit, the reconnection should be made by qualified service personnel.

■ Power Supply:	UL/CUL	TUV/UK
Volts:.....	110V AC	230V AC
Frequency:	60HZ	50 HZ

■ Product Specifications

Volts:.....	24V DC	24V DC
Amps:	750 mA	800 mA
Approximate rotational speed:	33 to 76 rpms	34 to 71 rpms
Magnet:.....	1,000 gauss	1,000 gauss
Magnet strength at tip of balls:.....	624 gauss	624 gauss
Magnet strength at rotator plate:	900 gauss	900 gauss
Polarity:	Negative field	Negative field

Benefits

The Magnassager combines the benefits of massage, magnetic therapy, and acupressure. The actual benefits are specific to the type of technique selected when applying the Magnassager. Because the Magnassager is ergonomically designed, its weight and rotating mechanism allows the user better control while performing the massage. By applying uniform strokes and constant pressure, the effectiveness of massage is greatly enhanced. In addition, forearm, hand and finger fatigue which is a common experience with hand massage is avoided when properly applying the Magnassager. Often, the root of muscle discomfort lies deep within the tissue layers and requires a progressive, thorough massage technique. This is one reason why the Magnassager is so effective.

When to Use and When NOT to Use the Magnassager

Indications

The Magnassager is designed to assist in providing a more thorough massage of the soft tissues, and is best indicated for:

- Tightness
- Stiffness
- Soreness
- Low back pain
- Fatigue
- Tension
- Muscle spasm
- Chronic back pain

The Magnassager is not recommended for individuals with the following conditions. If you have any of the following conditions consult with your physician before using this appliance.

Individuals with pacemakers and internal defibrillators should not use the Magnassager.

- Cancer
- Swelling that might be a malignancy
- Skin infections or ulcerations
- Edema or inflammation whether superficial or deep, trauma or hemorrhage (i.e., varicosities, peptic ulcers, hemophilia, aneurysm)
- Ischemia
- Bleeding disorders
- Phlebitis, thrombosis or other vascular diseases
- Lymphangitis
- Acute fever
- Pregnancy
- Acute osteomyelitis, tubercular joints or any infectious bone or joint disease
- Hypertension (do not perform vigorous massage on someone with hypertension)
- Intermittent claudication or symptoms of pain in the legs which gets worse upon exertion
- Insensitivity or any conditions where there is impairment of nerve sensation or an inability to communicate and/or respond to pain

Use of the Magnassager should be pleasant and comfortable. If the individual exhibits pain or discomfort, discontinue use immediately and consult a physician.

There are certain individuals in which the Magnassager should be used with caution and/or is contraindicated. Please consult a physician if you are unsure of its application.

Special Considerations

■ Pacemakers

The Magnassager should not be used on individuals with pacemakers and/or internal defibrillators.

■ Pain

The movement of oxygen-rich blood into the muscle, as well as the removal of toxins, may result in discomfort for some individuals. The user should be aware of the individual's tolerance and comfort level during the massage and adjust his/her technique accordingly.

■ Low Muscle Mass

Individuals with low lean body mass or those with loose or flaccid skin have bony protrusions and should use caution (e.g., scapula, iliac crest, ankles, elbows, and knees). Those persons with loose skin are at risk because the skin may collect under the massaging balls. The user should use their guiding hand to stretch the skin under the unit to avoid this situation.

■ Hair

Although unlikely, the Magnassager may pull or tangle body hair. This may be uncomfortable for the individual, therefore the user should encourage feedback or choose to perform the massage over an article of clothing. A water-based topical gel, lotion or oil over skin helps to inhibit pulling.

■ Alcohol

Alcohol impedes the massage process and may interfere with the elimination of toxins. Alcohol should be avoided within 24 hours of massage and directly after. Alcohol effects may mask pain sensors and may result in further injury or discomfort.

■ Diabetes

Diabetic individuals may have peripheral neuropathy and circulatory limitations. Massaging the feet and legs should be performed with caution.

■ Pregnancy

The Magnassager should not be used during pregnancy.

■ Medications

Certain prescription drugs (i.e., pain medication, heart medication, blood pressure medication, corticosteroids and blood thinners) may mask the pain sensors and require caution.

■ Dehydration

Because massage stimulates lactic acid, carbon dioxide, lymph, blood and oxygen flow, the user should encourage liberal water intake after each session.

How to Use the Magnassager

■ Breaking in period

The Magnassager motor may require a breaking-in period for up to 20 hours of use.

■ Preparation

Prepare for massage by creating a relaxing environment, encourage the individual to breathe during the treatment and, when possible, apply a thin layer of body gel, lotion or oil to the area to be massaged.

■ Technique

Note: Start with a slow speed to get used to the Magnassager's action. Stay in the muscle; let the Magnassager do the work for you.

With the Magnassager, the user should first massage the entire selected muscle group or muscle belly. Set Magnassager at a moderate to fast rotation and stroke the intended surface area at a speed of approximately two to three inches per second. The massage technique should include a general massage which outlines the muscle and the massage strokes should be longitudinal (i.e. up and down the length of the muscle.). During this technique, the user should let the weight and rolling action of the Magnassager glide over the skin. Avoid resting the unit in one place for an extended period of time. A gradual, progressive massage helps to relax the individual.

The initial technique should progress to a slower rotational speed and the massage strokes should attempt to cover a smaller surface area at a speed of approximately one inch per second. The massage should continue to cover the length of the muscle and allow the massage balls to transverse or massage across the muscle belly. Strokes should overlap and the user should attempt to massage the peripheral fluid in the direction toward the heart. This can be accomplished by applying firm pressure toward the heart and softer pressure away from the body. This procedure enhances the benefits received from the Magnassager.

Helpful Hints

- If the massage balls are cold, remove them from the unit, and place them under warm tap water for 5 seconds. NEVER RUN WATER OVER THE BALLS WHILE THEY ARE IN THE UNIT. The stainless steel massage balls are an excellent conductor of heat. Be sure to dry the balls thoroughly before placing them back into the unit.
- Better results can be obtained whenever the muscle can be elongated or stretched (e.g., ankle flexion to stretch the calf muscles). Use caution with recent or chronic injuries to the muscles(s) or skeleton.
- Always keep the Magnassager perpendicular to the body so that all three massage balls contact the skin evenly. Avoid bony protrusions; massage only the muscle.
- When using over clothing, use your free hand to keep the material taut to avoid binding.

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Copyright © 2002 by Thane International, Inc. All rights reserved.
Thane Fitness® is a division of Thane International, Inc.
Thane FitnessMD est une filiale de Thane International, Inc.
Droit d'auteur ©Thane International, Inc. 2002. Tous droits réservés.
Thane Fitness® es una división de Thane International, Inc.
Derechos de autor © 2002 por Thane International, Inc. Derechos reservados.